NEXT GENERATION VASCULAR ASSESSMENT SYSTEM

BASIC

VASCULAR DOPPLER PHOTO-PLETHYSMOGRAPHY







The newest technology in a space-saving design

Ideal for clinics or hospitals

CONFIGURATIONS


BASIC 1: Doppler

• 4 MHz & 8 MHz probes

BASIC 2: Doppler + bilateral photo-plethysmography

- 4 MHz & 8 MHz probes
- Two arterial photo-plethysmography sensors
- Two venous photo-plethysmography sensors





Doppler




Photo-plethysmography

Image: Introduction White background Image: Introduction ATTS MEDICAL WORLD Image: Introduction Image: Introduction

Atys medical



17, Parc d'Arbora 69510 SOUCIEU EN JARREST - FRANCE Tel. : +33 4 78 05 69 69 - Fax : +33 4 78 05 69 60 E-mail : atys@atysmedical.com http//www.atysmedical.com

BASIC Taking your Doppler device one step further

All the functions are fully integrated in the device. The BASIC is intuitive and easy to use.

DOPPLER

- 4 MHz & 8 MHz continuous Doppler
- Bidirectional Doppler
- Arterial & venous Doppler
- Real-time colour spectral analysis
- Audio output:
 - Two loud speakers
 - Stereo, separated reverse and forward flows
- Six indices displayed in real time (Vs, Vd, Vm, Hr, RI, PI, S/D)
- Display of the curve or spectrum

PHOTO-PLETHYSMOGRAPHY

- Ultrasensitive PPG sensors
- Record of the arterial waveform from the digits
- Diagnosis of venous insufficiency

USER INTERFACE

- 7" WVGA touch screen (TFT colour)
- Foot switch pedal

RESSOURCES

Buffer scrolling Embedded PC Storage of the examination reports Direct and high-quality printing of full-page colour reports Export of the examination reports to a USB memory stick Vessel labelling Patient data input

The BASIC is a portable device.

1000

The BASIC can be placed on an easily movable support called trolley I. This ergonomic trolley is fitted with a basket for storing accessories. Another trolley, named trolley Z, can support both the BASIC and an external printer.

> Dimensions (without trolley): 23 x 16 x 12 cm Weight (without trolley): 1.2 kg Power: 100V-240V, 50-60 Hz



